

# The Most Innovative Treatment Procedure Takes a New Turn



## Introducing Ancor Pro Orthodontic Anchorage System



Now the ability to enhance treatment capabilities and reduce chair time is in your hands. With the turn of a screw, the new *Ancor Pro* Anchorage System from Ortho Organizers helps you gain absolute anchorage in locations previously impossible without invasive procedures or extraoral appliances - providing greater control of tooth

movement, increased patient comfort and decreased reliance on patient compliance.

*Ancor Pro* anchors provide unmatched benefits - they are self tapping, allowing for easy placement and removal chairside. The multi-functional top of the anchor offers maximum versatility and performance, plus the anchors are available in multiple lengths. Additionally, the *Ancor Pro* Orthodontic Anchorage System is backed by the rigorous quality control you've come to expect from Ortho Organizers - *providing confidence at every turn.*

**\$1,450 + GST**

*Order now for delivery late July*



### Kit Components:

- 4 x Anchors 6mm
- 4 x Anchors 8mm
- 4 x Anchors 10mm
- 1 x Screwdriver
- 1 x Long driver insert
- 1 x Short driver insert
- 1 x Drill
- 1 x Latchlock driver
- 1 x Thumb wheel
- 12 x X-Ray placement guides
- 1 x Autoclavable storage container

**Free Call Aust 1800 645 813 : NZ 0800 440 883**

**www.orthoorganizers.com.au : sales@orthoorganizers.com.au**

**ORTHO<sup>2</sup>**  
Ortho Organizers  
Australia and New Zealand

# Control and confidence at every turn.

## Superior versatility, simplicity and ease of use – by design.

Developed by a team with more than 25 years of implant industry experience, the *Ancor Pro* Orthodontic Anchorage System has been specially designed to provide anchorage in a wide range of clinical applications, while offering simplicity, ease of use and patient comfort.

### Self-tapping & self-drilling mechanics

Sharp screw tips and threads allow for easy placement and removal chairside without general anesthesia.\*

### Multi-functional single head

Supragingival head features an upper button and undercut for the attachment of elastic chain, closed loop coil springs and other auxiliary orthodontic devices. Tapered conical top features a lumen to accept wires up to .022" in diameter. This innovative multi-functionality means greater procedural versatility and streamlined inventory control.

### Variety of popular sizes

Available in three lengths (6mm, 8mm, 10mm), and 1.6mm diameter; made from Grade 5 Titanium; providing maximum strength for ultimate performance.

\*Local anesthesia is suitable, in most cases



**Free Call Aust 1800 645 813 : NZ 0800 440 883**  
**www.orthoorganizers.com.au : sales@orthoorganizers.com.au**

**ORTHO<sup>2</sup>**<sup>TM</sup>  
Ortho Organizers  
Australia and New Zealand