



# Maestro™

## 1<sup>st</sup> Molar Buccal Tubes

Classic orthodontic design with a smooth modern flare.

The *Maestro* Buccal Tube is a non-convertible, single-tube that harmonizes classic orthodontic designs with the needs of today's practitioners and patients. This design provides the precise control needed throughout all phases of treatment for a predictable finish. The *Maestro* Buccal Tube's large funneled entrance, side grip areas and color-coding help simplify treatment by allowing for easier identification, placement, positioning and wire insertion. The ease of use associated with this product also reduces chair time which makes the *Maestro* Buccal Tube an extremely cost-effective choice for your practice.

A U.S. manufacturer, *Ortho Organizers*® is dedicated to providing world class products and solutions to meet your every orthodontic need. As a full service partner we offer the tools to help expand your practice and diversify your treatment options.



**Product characteristics:**

**Large funneled entrance** acts as a guide for easy wire insertion, resulting in reduced chair time

**Lowest possible profile** to optimize patient comfort without hindering the precision and control needed for case completion

**Easy Placement –**

**Color-coded ID mark** for easy identification

**Side grip areas** for tweezers to hold and position bracket

**Defined buccal indent** for accurate stability and placement

**Compound contoured base** design for superior fit and bond strength

**Patient comfort** is ensured with smooth rounded surfaces and lowered mesioingival hooks

**80 gauge mesh** bondable pad for consistent and proven bond adhesion

Prescriptions: Roth and McLaughlin, Bennett and Trevisi Rx's in both .018" and .022" sizes

Available as a direct bond and weldable

**Metal Injection Molded (MIM) Stainless Steel** for strength and proven performance

McLaughlin, Bennett & Trevisi	Description	Color-Code	Torque	Angulation	Offset	Catalog No.			
						.018		.022	
						Loose	Direct Bond	Loose	Direct Bond
STD	Upper right	Purple	-14	0	10	902-203	902-203DB	902-303	902-303DB
STD	Upper left	Purple	-14	0	10	902-204	902-204DB	902-304	902-304DB
STD	Lower right	White	-20	0	0	902-253	902-253DB	902-353	902-353DB
STD	Lower left	White	-20	0	0	902-254	902-254DB	902-354	902-354DB

Roth	Description	Color-Code	Torque	Angulation	Offset	Catalog No.			
						.018		.022	
						Loose	Direct Bond	Loose	Direct Bond
STD	Upper right	Purple	-14	0	14	902-201	902-201DB	902-301	902-301DB
STD	Upper left	Purple	-14	0	14	902-202	902-202DB	902-302	902-302DB
STD	Lower right	Red	-25	0	4	902-251	902-251DB	902-351	902-351DB
STD	Lower left	Red	-25	0	4	902-252	902-252DB	902-352	902-352DB

For more information on our products, orthodontic laboratory services and educational offerings, please contact us:

Customer Service: 888.851.0533

In the U.S.: 800.547.2000

Outside the U.S.: +(1) 760 448 8600

Canada: CERUM 800.661.9567

To fax an order: 800.888.7244

To fax an order outside the U.S.: +(1) 760 448 8607

To email an order:

USASales@OrthoOrganizers.com

To email an order outside the U.S.:

IntlOrders@OrthoOrganizers.com



1822 Aston Avenue  
Carlsbad, CA 92008  
USA

OrthoOrganizers.com