



NeoLucent™

Ceramic Bracket System

Performance you'll notice, brackets your patients won't.



Actual Patient

NeoLucent Ceramic Bracket System

The *NeoLucent* Low Profile Ceramic Bracket System delivers superb aesthetics by blending with the natural tooth tone for a virtually invisible appearance. A straight wire bracket system, *NeoLucent* Ceramic Brackets are extremely strong and provide the definitive performance you demand.

A US manufacturer, *Ortho Organizers*® is dedicated to providing world class products and solutions to meet your every orthodontic need. As a full-service partner, we offer the tools to help expand your practice and diversify treatment options. So when your patients desire a cosmetic alternative, we have the right choice for you, and them.

Patient Satisfaction

Invisible appearance – Crafted to the highest aesthetic standards

Translucent, natural appearance from 99.9% pure polycrystalline Alumina


Stain resistant throughout treatment

Comfort – Designed with your patient's comfort in mind

Low profile (see U1 profile measurements below)

Smooth, polished surface


Rounded corners



U1 Profile Measurements



NeoLucent Ceramic Bracket Profile



Clinical Performance

Accurate placement achieved with an easily removable color-coded ID mark and rhomboid shape

Stronger bonding and predictable debonding established with the *Combo-LOK™* Base which utilizes optimum mechanical and chemical retention

Easy ligation with ample tie wing clearance

Superior sliding mechanics via a precision cut archwire slot

Breakage resistant with unsurpassed strength

Antiallergenic for sensitive patients as bracket does not contain nickel, chromium or other metals

Routine protocol with no special instruments or adhesives required

Finishing Touches


To ensure excellence in every one of your cases, we offer an extensive line of aesthetic auxiliaries compatible with the *NeoLucent Ceramic Bracket System*.

Aesthetic Archwire Options



Tooth Tone Plastic-Coated Archwires

Stain and crack resistant this archwire features a tooth-colored coating that blends with natural dentition and reduces friction.



Micro-Coated Bio-Kinetix™ Thermally Activated Nitanium® Archwires

The labial side of this archwire has a tooth colored micro-coating to enhance aesthetics.



Pearl Tone® Preformed Ligation Wires


Available in a white/pearl hue, this archwire features a super smooth nonstick coating.

Aesthetic Elastomeric Options



Safe-T-Tie™ Ligatures

Ligatures for individual patient use, avoiding cross-contamination. We offer an array of color options that blend with the bracket and dentition. Available in clear, white, tooth tone, natural glow and smoke.




Chain Elastic

Made from super-elastic long lasting material, these ligatures deliver gentle, continuous force for more predictable tooth movement. Multiple colors to choose from: clear, white, tooth tone, natural glow and smoke.



Light Force Chain Elastic

Delivers a 20% lighter force than our regular chain making it an excellent choice for low force treatment mechanics. Available in clear.



"The response of my patients has been outstanding — they have been very impressed with the look and feel of the NeoLucent Ceramic Brackets."

R. Bruce McFarlane

D.M.D., B.S.c.D., M.C.I.D., F.R.C.D.(C)

Diplomate: American Board of Orthodontics

Orthodontic Instruments



O² Select™ Anterior Debonding and Scaler

Quality construction and design allows for safe and easy removal of direct bond brackets as well as the remaining adhesive material from the tooth surface without patient discomfort. The instrument's precision edges work for metal and aesthetic brackets. Catalog No. 203-104











O² Select Angulated Debonding Pliers

Superior strength tungsten carbide working edges allow for easy removal of brackets and excess bonding material with a conveniently angled head. This instrument is a great choice for removing anterior and posterior brackets. Catalog No. 203-135

NeoLucent Brackets









McLaughlin, Bennett and Trevisi Prescription*

.022		DESCRIPTION					
700-399MC		Set with hooks on 3s, 4s and 5s					
700-399MCQ		Set with hooks on 3s, 4s and 5s, single patient set					
SET (5X5)		COLOR-CODE	TOOTH NUMBER	TORQUE	ANGULATION	ITEM NUMBERS	
MAXILLARY						.022 R	.022 L
Centrals		black	U1	+17°	+4°	700-307	700-308
Laterals		pink	U2	+10°	+8°	700-323	700-324
Cuspids w/Hook		green	U3	-7°	+8°	700-339HK	700-340HK
1 st & 2 nd Bicuspids w/Hook		purple	U4&5	-7°	0°	700-357HK	700-358HK
MANDIBULAR						.022 R	.022 L
Anteriors		yellow	L1&2	-6°	0°	700-331	700-332
Cuspids w/Hook		blue	L3	-6°	+3°	700-371HK	700-372HK
First Bicuspids w/Hook		gray	L4	-12°	0°	700-385HK	700-386HK
Second Bicuspids w/Hook		white	L5	-17°	0°	700-387HK	700-388HK

Please note all color-code dots are gingival distal (for brackets with a single dot).

NeoLucent Brackets

Roth Prescription*

.018		.022		DESCRIPTION						
		700-397	Set with no hooks							
		700-397Q	Set with no hooks, single patient set							
700-298	700-398	Set with hooks on 3s								
700-298Q	700-398Q	Set with hooks on 3s, single patient set								
700-299	700-399	Set with hooks on 3s, 4s and 5s								
700-299Q	700-399Q	Set with hooks on 3s, 4s and 5s, single patient set								
SETS (5X5)		COLOR-CODE	TOOTH NUMBER	TORQUE	ANGULATION	ITEM NUMBERS				
MAXILLARY						.018		.022		
						.018 R	.018 L	.022 R	.022 L	
Centrals		black	U1	+12°	+5°	700-205	700-206	700-305	700-306	
Laterals		pink	U2	+8°	+9°	700-221	700-222	700-321	700-322	
Cuspids		green	U3	-2° -2°	+11° +11°	700-235 700-235HK	700-236 700-236HK	700-335 700-335HK	700-336 700-336HK	
1st & 2nd Bicuspid		purple	U4&5	-7° -7°	0° 0°	700-257 700-257HK	700-258 700-258HK	700-357 700-357HK	700-358 700-358HK	
MANDIBULAR						.018 R	.018 L	.022 R	.022 L	
Anteriors		yellow	L1&2	0°	0°	700-233	700-234	700-333	700-334	
Cuspids		blue	L3	-11° -11°	+5° +5°	700-267 700-267HK	700-268 700-268HK	700-367 700-367HK	700-368 700-368HK	
1st Bicuspid		white	L4	-17° -17°	0° 0°	700-287 700-287HK	700-288 700-288HK	700-387 700-387HK	700-388 700-388HK	
2nd Bicuspid		red	L5	-22° -22°	0° 0°	700-291 700-291HK	700-292 700-292HK	700-391 700-391HK	700-392 700-392HK	

Please note all color-code dots are gingival distal (for brackets with a single dot).

For more information on our products, orthodontic laboratory services and educational offerings, please contact us:
Customer Service: 888.851.0533
In the U.S.: 800.547.2000
Outside the U.S.: +(1) 760 448 8600
Canada: CERUM 800.661.9567
To fax an order: 800.888.7244
To fax an order outside the U.S.: +(1) 760 448 8607
To email an order: USASales@OrthoOrganizers.com
To email an order outside the U.S.: IntlOrders@OrthoOrganizers.com



1822 Aston Avenue
Carlsbad, CA 92008
USA
OrthoOrganizers.com